

AEROSPACE / INDUSTRIAL BIOPHARMACEUTICAL BIOTECHNOLOGY MEDICAL DEVICE PHARMACEUTICAL

FACE MASK Diamond Shape



Aerospace / Life Sciences / Microelectronics



FACE MASK Diamond Shape

Features

Designed to have more room for easier breathing by eliminating discomforts to the nose in exhaling.

Keep the face mask at a distance from the lips and mouth, providing more room for easier speech.

Face mask can be adjusted to any facial shape and inner space can be enlarged for greater comfort.





FACE MASK Diamond Shape

Applications

Specification

Suitable for emergency rooms, operating rooms, dental procedures, clinics, cleanrooms, laboratories, pharmaceutical industry, among others.

- Type: Diamond shape, 3-ply, non-sterile.
- Outer Layer: Water-repellent non-woven.
- Middle Layer: High level static Filtration nonwoven.
- Inner Layer: Micro-fibered Thermal bond nonwoven.
- · Nose band: PE Iron whire.

INDUSTRIAL SOLUTIONS

FACE MASK Diamond Shape

Caution

- Do not pull the body of the facemask to over from between the face and the mask.
- Replace with a new set of face mask upon entering or exiting potentially infectious area, and at times when either one become unusable (i.e. has a rip or a tear).
- Replace when there is unusual smell in breathing.
- Avoid contact with the outer surface of the face mask when taking it off to prevent infection and contamination.
- Store in a cool and dry place.





FACE MASK Diamond Shape

Contact information

Quotes and orders from around the world

APPS INC. appsindustrial.com +1 (818) 930 0779 [Toll free]: +1 818-5743-345 contact@appsindustrial.com

APPS INC.



Aerospace / Life Sciences / Microelectronics