



AEROSPACE / INDUSTRIAL  
BIOPHARMACEUTICAL  
BIOTECHNOLOGY  
MEDICAL DEVICE  
PHARMACEUTICAL

# ANTISTATIC FABRIC COVERALLS MEASUREMENT GUIDE

Tyvek and / or polypropylene  
(collared garments)

Life Sciences / Industrial



## ANTISTATIC FABRIC COVERALLS MEASUREMENT GUIDE

Tyvek and / or polypropylene  
(collared garments)

SIZE	CHEST CONTOUR (CM)	CHEST CONTOUR (INCH)	TOTAL LENGTH (CM)	LONG SLEEVE FROM THE CENTER OF THE BACK (INCH)	GUSSET (INCH)
36 CH	112.5	45	1.54	34.75	30
38 MED	117.5	47	1.56	36	32
40 LARGE	122.5	49	1.58	37.25	34
42 XLARGE	127.5	51	1.6	38.75	35
44 XXLARGE	132.5	53	1.62	40.25	36
46 XXXLARGE	137.5	55	1.64	41.25	37
48 XXXXLARGE	142.5	57	1.66	42.75	38
50 XXXXXLARGE	147.5	59	1.68	43	38

The slack of the garments is 9 inches, it is considered in each size, this is already taken into account in the contour of the thorax. Measurements are only a reference point.

ANTISTATIC FABRIC  
COVERALL  
MEASUREMENT GUIDE

Tyvek and / or polypropylene  
(collared garments)

# Contact information

Quotes and orders from around the world

APPS INC.  
appsindustrial.com  
+1 (818) 930 0779  
[Toll free]: +1 818-5743-345  
contact@appsindustrial.com

